Application and safety instructions for the MEISINGER Tenting Screw System

CONTENT

Language	Page
English (EN)	2-4
French (FR)	5-7
German (DE)	8-10
Italian (IT)	11-13
Spanish (ES)	14-16

Application and safety instructions for MEISINGER Tenting Screw System

Screwssystem

Art.-No.: 79BTS00

The Bone Management® Tenting Screw System is a treatment option for both simple and complex bone defects and offers an alternative to the shell technique and the use of titanium-reinforced membranes. The screws contained in the system are made of the titanium alloy Grade 5-ELI (Titanium Grade 23) and have a self-tapping thread. The screw head includes an interface to accommodate the insertion tool (hexalobular socket).

INDICATIONS FOR USE

The Tenting bone screws are indicated for the fixation and stabilisation of bone grafts and bone substitutes as well as for the support of barrier membranes. They are used in dental surgery therapies (osteosynthesis, correction of degenerative diseases) to restore anatomical structures, but are not suitable as a replacement for such structures.

UTILIZATION

The bone screws can be used for a variety of surgical techniques depending on the case situation. The selection of the technique and the composition of the instruments to be used must be selected and prepared for use before the operation begins. The user must ensure that the screws are secured against aspiration both during fixation and removal (e.g. by placing a rubber dam).

Preparation for screws

Due to the self-tapping thread design, no pilot hole is usually required for inserting the screws. If required, the bone can be pre-punched using the initial bur included in the system. Pre-graining is recommended for hard bone. In the case of very hard bone, pre-drilling using the drill included in the system may be indicated.

Selecting the screwdriver and holding the screws

Depending on the selected procedure (mechanical or manual screw insertion), the appropriate screwdriver from the manufacturer must be selected. The TL0T1 screwdriver is used for manual insertion. The 31053 screwdriver is used for automatic insertion with the contra-angle handpiece. The use of third-party accessories may impair the function of the screws and invalidates any warranty or replacement.

The screws are contained in a sterilisable box that can be opened using a sliding mechanism. The screws are removed from the box using the screwdriver. To do this, the screwdriver is inserted vertically into the interface contained in the screw heads. The screws are not fixed in the box. Therefore, care must be taken to ensure vibration-free working both when opening the box and when picking up the screws with the screwdriver. As the screws are implantable products, the screws must only be picked up using the screwdriver to avoid touching them with your hand.

Postoperative procedure

After clinically and radiographically confirmed, complete bone healing and functional restoration, the screws are either removed or left in situ, depending on the case situation and the aim of the procedure. If an infection or weakening of the bone tissue occurs in the area of the inserted bone screw, the screw must be removed. Explantation should also be performed in the event of screw loosening and always in children to avoid growth inhibition.

Sterilisation of the screws

The screws are supplied non-sterile and must be reprocessed before use (cleaning / disinfection / sterilisation). For proper reprocessing of the items, please observe the instructions for reprocessing (cleaning, disinfection and sterilisation) of medical devices from Hager & Meisinger GmbH.

Disposal

When disposing of the instruments (at the end of their service life or after expiry of the specified shelf life), ensure that the product is disposed of in the waste for biohazardous substances. All packaging components are disposed of in accordance with national regulations (e.g. dual waste system).

Contraindications

Basically, general medical as well as local, absolute and relative contraindications for dental surgical measures must be observed.

Absolute contraindications

- Active infections in or near the area to be augmented as well as local or systemic pathological processes (e.g., symptoms such as fever, local inflammation, abscesses)
- Insufficient bone supply (quality / quantity)
- Illnesses or medication that impair bone metabolism
- Insufficient soft tissue supply, which does not allow complete coverage of the augmented area

Relative contraindications

- Dentoalveolar growth not completed (exception: cases where no dentoalveolar growth is expected, e.g., ectodermal dysplasia) cases where no dentoalveolar growth is expected, e.g., ectodermal dysplasia)
- Known allergies and / or foreign body reaction to the implant material
- Systemic and/or metabolic illnesses or medical treatments that lead to progressive bone deterioration (e.g., cortisone, immunosuppressants, bisphosphonates and diabetes)
- Drug and alcohol abuse
- Lack of patient compliance
- Poor circulation
- Heavy work or active sports
- Mental state which may lead to non-compliance with the medical order
- Highly atrophied jaw
- Nicotine consumption, which may jeopardize the success of the procedure by delaying bone wound healing

IMPORTANT: Care must be taken to protect the anatomical structures (safety distance min. 2 mm) as well as the course of the adjacent teeth / tooth roots (hazard of damage, infection / resorption). To avoid the development of infections and dehiscence, care should be taken to prevent contact between the graft and adjacent teeth. Cases in which the mechanical stability and mechanical rest of the graft cannot be achieved may compromise successful bone healing.

SEPERATE INSTRUCTIONS FOR SCREWS: When selecting the screw length as well as the screw placement, the evaluation of the available bone supply by means of X-ray diagnostics must be taken into account. To prophylactically counteract soft tissue perforation or dehiscence, tension-free and complete wound closure must be achieved. The entire augmentation area, including the bone screws, must not be loaded during the healing process. Failure to do so may result in bending, loosening, loss and/or breakage of the screws and this means a negative impact on the augmentation result. Also, excessive physical activity as well as trauma affecting the augmentation site can lead to premature failure of the screws by loosening or breaking. The screws must only be used in accordance with the instructions for use and safety. The complete sinking of the screw can be limited by very hard bone quality. Any kind of modification in the procedure or on the product can lead to impairment of the treatment success or damage to the product. The screws are a single-use product and must therefore not be reprocessed and reused after use. Even screws that have already been used but are visually undamaged may suffer microscopic damage from previous mechanical stress leading to fatigue and product failure. The screws have not been evaluated for safety and compatibility in the MR environment. They were not tested for heating, migration, or image artefacts in the MR environment. The safety of the screws is uncertain in the MR environment. Scanning a patient with this medical device could result in injury to the patient.

Any serious-incident that has occurred in relation to the device should be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

The "Summary on Safety and Clinical Performance Report" (SSCP) is a summary of the clinical data and further information on the safety and clinical performance of the screws and can be provided as a copy upon request by e-mail to RA@meisinger.de. Once the European Database for Medical Devices (EUDAMED) is operational, the SSCP for the screws can be viewed on the following website: <https://ec.europa.eu/tools/eudamed>.

General instructions

Please follow general application and safety instructions for MEISINGER products in the medical area and also the advice for processing (cleaning, disinfection and sterilisation) of medical devices from Hager & Meisinger GmbH.

Manufacturer

Item number

Medical device

Batch number

Do not re-use

Non-sterile

Under U.S. federal law, this device may only be sold to or on behalf of trained medical professionals.

Caution

Conformity mark

Do not use if the packaging is damaged

Manufacturing date

Packaging unit

Consult electronic instructions for use


ifu.meisinger.de


Hager & Meisinger GmbH
 Hansemannstr. 10
 41468 Neuss • Germany
 Tel. +49 2131 2012-0
 Fax: +49 2131 2012-222
 Email: info@meisinger.de
 Internet: www.meisinger.de

CONTENT

	Screws		
			
	7	7	7
Fig.	TS035 040 030	TS035 040 040	TS035 040 050
L ₁ (mm)	3.5	3.5	3.5
L ₂ (mm)	4.0	4.0	4.0
L ₃ (mm)	3.0	4.0	5.0
L ₄ (mm)	7.75	8.75	9.75
	CE 0044 / 		

	Screws					
						
	3	3	3	2	2	2
Fig.	US050 040 070	US050 040 090	US050 040 110	US070 040 070	US070 040 090	US070 040 110
L ₁ (mm)	5.0	5.0	5.0	7.0	7.0	7.0
L ₂ (mm)	4.0	4.0	4.0	4.0	4.0	4.0
L ₃ (mm)	2.25	4.25	6.25	2.25	4.25	6.25
L ₄ (mm)	7.0	9.0	11.0	7.0	9.0	11.0
	CE 0044 / 					

	Initial bur	Twist drills	Screwdriver TX RA long
			
	1	1	1
Fig.	186RF	203RF	31053
Schank ¹	204	205	205
Size ²	018	011	-
Length mm	12.0	9.0	27.0
	-	-	-
	1.8	1.1	-
Opt. speed rpm	800	1000	10
Max. speed rpm	6000	6000	15
	CE 0044		

	Screwdriver
	
	1
Fig.	TL0T1
	Figure 1:2 (length mm: 148.0)
	CE

¹ 204=RA, 205=RA L ² Largest working part diameter in 1/10 mm  Minimal diameter  External diameter

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