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Implantate A = Rot = Ø 3,5mm
Implantate M = Orange = Ø 4,0mm
Implantate B = Gelb = Ø 4,5mm

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Preparation of the instruments


The instruments and the surgical cassette are not supplied sterile. Before each usage, the instruments must be checked for their suitability for use. Please observe that the "General Instructions for Use and Safety for MEISINGER products in the medical field and also the instructions for reprocessing (cleaning, disinfection and sterilisation) of medical devices from Hager & Meisinger GmbH". Care must be taken to use the instruments specifically suited to the implant variant.

Apical overlengths

When preparing the implant bed, it must be taken into account that the effective drilling depth varies by 0.4 - 0.7mm apically and, depending on the diameter of the implant drill, is deeper than the targeted implant length. This additional length must already be taken into account during the planning phase.

	Tri-Spade Drill A	Tri-Spade Drill M	Tri-Spade Drill B
Diameter	Ø2.4mm	Ø2.9mm	Ø3.3mm
Apical overlengths	0.4mm	0.5mm	0.6mm


Apical overlength using the example:


Explanation of the laser marking

The laser markings on the twist drills and Tri-Spade drills are based on the nominal dimensions of the implant lengths. The lower and the upper edge of a marked block correspond to one measurement each. Starting with the tool tip, they match the implant lengths 6.6/8,0/9,5/11/14.


Using the Tri-Spade drill as an example:


Using the example of Conical Reamers:

Using the example of Cortical Countersink:

Using the example of a tap:

Comment: Use a depth gauge with legal SI metric units to determine the exact depth of the implant bed. The laser markings of the drills/conical reamers/cortical countersinks/taps do not serve to accurately measure the drilling depth.

Recommended speeds/torques Instrument set Surgery Cassette

Instrument	Recommended number of revolutions	Recommended torque
Rose drill	2.000rpm	-
Initial drill	1.000rpm	-
Twist drill	800rpm	-
Tri-Spade Drill	800rpm	-
Conical Reamer	50rpm	max. 50Ncm
Cortical Countersink	50rpm	max. 50Ncm
Tap	15rpm	max. 50Ncm
Myplant implant	15rpm	max. 50Ncm

Surgical procedure Instrumentation Surgery Cassette - Surgical procedure with myplant implants

Implant	Rose drill	Initial drill	Twist drill	Tri-Spade Drill A Ø2.4	Tri-Spade Drill M Ø2.9	Tri-Spade Drill B Ø3.3	Tri-Spade Drill M Ø2.9	Tri-Spade Drill B Ø3.3	Tri-Spade Drill M Ø2.9	Tri-Spade Drill B Ø3.3	Conical Reamer A (right-handed/cutting)	Conical Reamer B (counter-clockwise/condensing)	Conical Reamer M (counterclockwise/condensing)	Conical Reamer B (clockwise/cutting)	Conical Reamer M (clockwise/cutting)	Conical Reamer B (counterclockwise/condensing)	Cortical Countersink A	Cortical Countersink M	Cortical Countersink B	Tap A	Tap B
Implant A Ø3.5mm	x	x	x	x	(X)						(X)	(X)									
Implant M Ø4.0mm	x	x	x	x	x	(X)					(X)	(X)									(X)
Implant B Ø4.5mm	x	x	x	x	x	x	(X)				(X)	(X)									(X)

X = always recommended
(X) = always recommended for harder bone quality (e.g. D1/D2)
(X) = always recommended for softer bone quality (e.g. D3/D4)
O = Optional use after Tri-Spade drilling for very soft cancellous /bone quality. Conical Reamer and Cortical Countersink do not have to be used in this case.

Comment:

- For gentle reprocessing, the instruments should be used according to the recommended number of revolutions and under low pressure. Ensure that the instruments are sufficiently cooled during preparation.
- Any contact and contamination with the implant must be avoided.
- It should be noted that the twist drills have an apical overlength at the drill tip, i.e. the preparation depth of the implant bed does not correspond to the countersinking depth of the implant. This must be taken into account when choosing the implant length.
- Patient evaluation/preoperative treatment planning: It is essential that these steps precede the surgical part, as inadequate planning can lead to implant loss.
- Opening/exposure: The gingiva is opened with a scalpel/gingiva punch. This can be followed by the preparation of a mucoperiosteal flap to visualise the osseous structures.
- Rose drill: For smoothing the alveolar ridge, removing unevenness
- Initial drill: Initial drilling/marking hole
- Twist drill: Alignment and first deep drilling according to the implant length.
- Directionality: Used to check the implant axis and parallelism if necessary. To check the alignment of the first deep borehole, the slim end is placed in the borehole.
- Tri-Spade Drill: For preparing the implant bed according to implant diameter and bone quality. Tri-Spade drills in two diameters are available for each implant diameter (A/M/B). If the drilling tool with the smaller diameter (EPS, marked with a colour ring for colour coding) is used for the final Tri-Spade drilling, the aim is to achieve increased primary stability. This procedure is recommended for very soft bone qualities (e.g. D3/D4). When using the more diameter Tri-Spade Drill (RPS, marked with two colour rings for colour coding) for final Tri-Spade drilling, regular primary stability is aimed for. This procedure is recommended for harder bone qualities (e.g. D1/D2).
- Directionality: Used to check the implant axis and parallelism if necessary. To check the alignment and depth of the implant bed after Tri-Spade drilling, the wide end is placed in the drill hole.
- Conical Reamer in clockwise rotation: For conical preparation of the implant bed with harder bone qualities (e.g. D1, D2).
- Conical Reamer in counter-clockwise rotation: For conical preparation of the implant bed with softer bone qualities (e.g. D3, D4) to compress/condense the bone.
- Cortical Countersink: Optional, for widening the cortical bone, in case of very soft cancellous bone structure, analogous to the diameter of the implant.
- Tap: For pre-cutting the implant thread. If the bone quality is compromised (e.g. D3, D4), the tap can be dispensed with.

Open-end spanner/ Sleeve spanner

After inserting an implant, the insertion unit is removed. To remove them, the open-end spanner (open-end spanner 3,0 (small)) is placed over the narrowest part of the insertion instrument on the hexagon of the insertion unit and locked. To loosen the screw of the insertion unit with the insertion instrument, turn the insertion instrument counter-clockwise. If the hexagon is not accessible (e.g. in narrow gaps), the sleeve spanner is used. It is placed on the mechanical insertion instrument and then on the hexagon of the insertion unit before the insertion process. The sleeve spanner has a hexagon at the upper end. The larger jaw of the open-end spanner (open-end spanner 3,5 (large)) is now used for counter-holding. The sleeve spanner thus moves the point for counter-holding (loosening the insertion unit) above the level of the occlusal surface of the adjacent teeth and thus makes the hexagon easier to reach.

Further information can be found in the myplant manual (85MMP2) and at www.myplant-dental.com for further information.

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Implantes A = rojo = Ø 3,5mm
Implantes M = naranja = Ø 4,0mm
Implantes B = amarillo = Ø 4,5mm

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