DE

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EN

Automatic, Cleaning and disinfection device:

Securely apply the cleaned ratchet parts to the carrier. Do not overload the carrier. Start the program. After rinsing, chemical cleaning starts at 40°-60°C. Residues from the cleaning process must be reliably removed in the subsequent rinsing phase. Avoid damage to the material from neutralising reagents. Thermal disinfection is achieved at 90°-95°C. The subsequent treatment with deionised water is followed by adequate drying. Remove the ratchet parts from the device immediately after the program ends.

Maintenance, inspection and testing

Allow the parts to cool to room temperature and visually inspect them for residues of proteins and other soiling. If necessary, repeat the processing steps.

Lightly lubricate the areas marked with ⚠ using contra-angle handpiece oil. Assemble ratchet and carry out functional test.

Sterilisation packaging

Place the ratchet into packaging suitable for sterilisation according to ISO 11607 and EN 868. The bag must be large enough for the instrument. Closures must not be under tension.

Sterilisation

Method: Fractionated pre-vacuum process in accordance with ISO 17665

Temperature: heating to 134°C, max. 137°C

Pressure: 3 pre-vacuum phases with a minimum pressure of 60 millibars

Holding period: 5 minutes

Drying time: at least 10 minutes

After sterilisation, inspect the sterile packaging for damage and check sterilisation indicators.

Place of use

Immediately after use, the products must be put in cold water (<40 °C). Do not use hot water (> 40 °C) or cleaning agents, as this may cause fixation of residues on the product (risk of protein coagulation (denaturation)), which can affect the success of the subsequent processing steps.

Storage

Store the ratchet at a moderate temperature, and in a dry, dust-free, well-ventilated place, in which there is no corrosive steam.

Adjustment / check

We recommend annual adjustment / check of the ratchet. The ratchet must be cleaned and sterilised before shipping, otherwise the product will be returned.

General instructions Please follow general application and safety instructions for MEISINGER products in the medical area and also the advice for processing (cleaning, disinfection and sterilisation) of medical devices from Hager & Meisinger GmbH. Further information The preparer is responsible for ensuring the desired details: using the equipment, materials and personnel. The process must be validated and routine monitoring of the process is required with used equipment. Further information on the preparation of medical devices at www.rki.de or www.a-k-i.org ATTENTION: Any serious-incident that has occurred in relation to the device should be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.	
---	--

									
Manufacturer	Sales partner	Item number	Medical device	Batch number	Non-sterile	Follow the instruction for use	Under U.S. federal law, this device may only be sold to or on behalf of trained medical professionals.	Caution, consult accompanying documents	Manufacturing date

	Hager & Meisinger GmbH Hansenannstr. 10 41468 Neuss • Germany Tel. +49 2131 2012-0 Fax: +49 2131 2012-222 Email: info@meisinger.de Internet: www.meisinger.de
	Meisinger USA, L.L.C. 10150 E. Easter Avenue Centennial, Colorado 80112 • USA Phone: +1 (303) 268-5400 Fax: +1 (303) 268-5407 Email: info@meisingerusa.com Internet: www.meisingerusa.com

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