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Abkürzung	Bedeutung
TL RP	Tissue Level Regular Platform
TL WP	Tissue Level Wide Platform
BL NC	Bone Level Narrow Connect
BL RC	Bone Level Regular Connect

Application and safety instructions for MEISINGER MEISINGER

OKTAGON Prosthetic components

These instructions for use apply to the cover screws, gingiva formers, cementable abutments, temporary abutments, solid abutments, screw-retained abutments, titanium adhesive bases, ALLIGATOR abutments, ball attachments, pre-milled titanium blocks and Multi-Unit abutments in the OKTAGON implant line.

INTENDED USE

The prosthetic components are indicated for use in the restoration of dental implants during the healing phase after implantation, for enabling subgingival and transgingival healing, and for enabling prosthetic restoration in the upper and lower jaw for various indications on the basis of a dental implant.

SCOPE OF APPLICATION

The OKTAGON prosthetic components are suitable for the restoration of OKTAGON dental implants of different types, such as subgingival and transgingival implants with different lengths, diameters and platforms. Prosthetic applications may include individual crowns, bridges, partial or full dentures. Depending on the indication and anatomical conditions, the user selects the various components for the individual prosthetic purpose.

GENERAL PRODUCT DESCRIPTION

The choice of prosthetic restoration and the prosthetic components required for it has a significant influence on the long-term success of the implantation. Intensive communication between the dentist and dental technician, careful pre-prosthetic planning and the involvement of the patient are important pre-requisites for this. Prosthetic treatment can be carried out using single crowns, bridges, partial or full dentures and can be cemented, screwed or integrated as a removable prosthesis on the corresponding abutments for OKTAGON implants.

The OKTAGON implant line offers a wide range of different prosthetic components for the patient's individual prosthetic treatment. These include cover screws, gingiva formers, cementable and screw-retained abutments, solid abutments, titanium adhesive bases, ALLIGATOR abutments, ball attachments, pre-milled titanium blocks and Multi-Unit abutments. The components are connected to the implant using a metric thread and, in some cases, a rotation-locked interface.

PRODUCT-SPECIFIC INSTRUCTIONS FOR THE PROSTHETIC PROCEDURE

The following table shows which prosthetic components are available with or without rotation locking.

Prosthetic component	Design with rotation-locked interface	Design without rotation-locked interface
Cover screws		✓
Gingiva formers		✓
Cementable and screw-retained abutments	✓	
Temporary abutments	✓	
Titanium adhesive bases	✓	✓
Ball attachments		✓
Solid abutments		✓
Pre-milled titanium blocks	✓	
Angled ALLIGATOR abutments	✓	
Straight ALLIGATOR abutments		✓
Angled Multi-Unit abutments	✓	
Straight Multi-Unit abutments		✓

All prosthetic components and instruments are supplied non-sterile and are labelled with the appropriate symbols. They must be cleaned, disinfected and sterilised before use.

All prosthetic components are designed for single-use. Reusing single-use products carries an increased risk of infection and risk-free functional safety (e.g. in terms of stability or accuracy of fit) cannot be guaranteed.

COMPATIBILITY INFORMATION

OKTAGON implants must always be used with OKTAGON prosthetic components. When doing so, it is important to ensure compatibility between the implant and the prosthetic component by only selecting components that match the implants to be restored (TL RP, TL WP, BL NC, BL RC). In addition, care must be taken to ensure that the prosthetic components are only used with the insertion aid provided.

COVER SCREWS AND GINGIVA FORMERS

After insertion, the OKTAGON implant is sealed. For subgingival healing, a sterile cover screw is included with each implant (TL: gingival height 0.0 mm, BL gingival height 0.5 mm). Alternatively, a cover screw with a different gingival height can be used, which is available separately. In the case of subgingival healing, a second procedure is necessary. After the healing phase, the gingiva is surgically opened at the implant site, the cover screw is removed and replaced with a suitable abutment. For transgingival healing, a gingiva former with the desired gingiva height is to be selected. To shape the soft tissue, the wound edges are sutured atraumatically close to the restoration. The suture closure must not impair the vascular supply of the mucosal edges. The selection of the appropriate cover screw or gingiva former is patient-specific and is the responsibility of the user. It is recommended to tighten the cover screw or gingiva former hand-tight (5-7 Ncm).

TEMPORARY IMMEDIATE RESTORATION

Temporary abutments can be used prior to placement of the final prosthetic restoration to maintain, stabilise and shape the soft tissue during the healing phase. It is recommended that the temporary abutment be tightened in the implant to a maximum torque of 35 Ncm.

A temporary prosthetic immediate restoration should only be carried out if the patient fulfils the requirements for osseointegration in the early phase; otherwise, a delayed solution is preferable. The temporary restoration should not be brought into occlusion.

IMPRESSION

The components used for the impression are hand-tightened (5-7 Ncm) in the implant. When using an impression material, care should be taken to select a suitable material. This is to be used according to the manufacturer's instructions.

The impression can be taken by means of open or closed impression-taking at implant level, by means of an impression cap at abutment level or digitally, so that a model with laboratory analogues can be produced afterwards. When using angled prosthetic components, it must be ensured that they are repositioned identically on the model in the patient's mouth (e.g. by observing the indexing of the prosthetic components or by fabricating a transfer key). When taking abutment-level impressions, non-rotationally secured prosthetic components must not be removed again after taking the impression.

PROSTHETIC TREATMENT

The following table shows the prosthetic treatment options.

Prosthetic component	Single tooth	Bridge	Partial denture	Total denture
Abutment	✓	✓	✓	✓
Solid abutment	✓	✓		
Titanium adhesive base	✓	✓	✓	✓
Ball attachment			✓	✓
ALLIGATOR			✓	✓
Pre-milled titanium block	✓	✓		
Multi-Unit abutment		✓	✓	✓

The final prosthetic restoration can be inserted and placed in occlusion as soon as the implant has sufficient stability or is fully osseointegrated. Angled abutments should not be used in regions with high mechanical loads on implants with a small diameter.

After the final restoration has been received, the previous temporary restoration is removed. The cleaned and sterilised restoration is then placed in the patient's mouth. Before tightening the screw, ensure that the implant-secondary part connection and the alignment are correct. Important: Attention should always be paid to the load distribution of the prosthetics.


Torques greater than 35 Ncm may cause abutment and/or implant failure. Torques below the recommended values may cause abutment loosening and subsequent abutment and/or implant failure. Once the screw has been tightened to 35 Ncm, it is recommended not to remove it again to prevent the screw from wearing out. If it becomes necessary to remove the restoration, it can be reattached using a new screw that matches the prosthetics.

The appropriate abutment is selected according to the implant used (TL RP, TL WP, BL NC, BL RC) and the planned restoration.

Abutment

In the case of angled abutments, in addition to the correct selection of the angle, the alignment of the angle must be observed. In the OKTAGON system, abutments are available in two versions – A and B. Version A is angled over the corner of the Oktagon, which is marked with a dot. Version B is angled over the surface of the Oktagon, which is marked with two dots.

The abutment is inserted using the corresponding system-compatible OKTAGON screwdriver. The abutment can be placed in 4 (BL) or 8 (TL) different positions. Care must be taken to ensure that the abutment is correctly seated in the implant interface. Taking into account the recommended torque of 35 Ncm, the abutment is now fixed in the implant by tightening the base screw using the screwdriver. It is recommended to check the abutment screw in the implant again after about a week, using the specified torque, and to tighten it if necessary.

Solid abutment

The desired height of the solid abutment is selected according to the planned restoration. The solid abutment is inserted using the corresponding system-compatible OKTAGON insertion tool. Care must be taken to ensure that it is correctly seated in the implant. The solid abutment is now tightened in the implant using the insertion tool, taking into account the recommended torque of 35 Ncm. It is recommended to check the specified torque again after about a week and to tighten it again if necessary.

Titanium adhesive base

The titanium adhesive base is inserted using the corresponding system-compatible OKTAGON screwdriver. The titanium adhesive bases with rotation lock can be placed in 4 (BL) or 8 (TL) different positions. Titanium adhesive bases without rotation lock can be placed freely in the implant. In all cases, make sure that the titanium adhesive base is correctly seated in the implant interface. Taking into account the recommended torque of 35 Ncm, the titanium adhesive base is now fixed in the implant by tightening the base screw using the screwdriver. It is recommended to check the screw in the implant again after about a week, using the specified torque, and to tighten it if necessary.

Ball attachment

The ball attachment is inserted using the appropriate system-compatible OKTAGON insertion tool. Care should be taken to ensure that it is correctly seated in the implant. The ball attachment is now tightened in the implant using the insertion tool, taking into account the recommended torque of 35 Ncm. It is recommended to check the specified torque again after about a week and to tighten it again if necessary.

ALLIGATOR

The appropriate ALLIGATOR is selected based on the implant used and the planned prosthetic restoration. The angled ALLIGATORS with rotation lock can be placed in 4 (BL) or 8 (TL) different positions. The orientation of the angle must be taken into account. The straight ALLIGATORS without rotation lock can be placed freely in the implant. In all cases, it must be ensured that the ALLIGATOR is correctly seated in the implant interface.

The ALLIGATORS are tightened using the following OKTAGON instruments:

- Straight: insertion aid (REF 75836)
- Angled: Torx screwdriver (REF 63978, 64551, 64553, 77900)

The recommended torque of 35 Ncm must be observed. It is recommended to check the specified torque again after about a week and to tighten it if necessary.

Pre-milled titanium block

The pre-milled titanium blocks are used to fabricate patient-specific CAD/CAM-milled abutments. All pre-milled titanium blocks are compatible with the Medentika holder and have a rotation-locked implant interface. The selection of the appropriate titanium block is made according to the diameter of the prosthetics to be fabricated. The manufactured abutment is inserted using the corresponding system-compatible OKTAGON screwdriver. The abutment can be placed in 4 (BL) or 8 (TL) different positions. The alignment of the angulation must be observed. Care must be taken to ensure that the abutment is correctly seated in the implant interface. Taking into account the recommended torque of 35 Ncm, the abutment is now fixed in the implant by tightening the base screw using the screwdriver. It is recommended to check the abutment screw in the implant again after about a week, using the specified torque, and to tighten it if necessary.

Multi-Unit abutment

The appropriate Multi-Unit abutment is selected based on the implant used and the planned restoration. The Multi-Unit abutments are placed in the implant using the pre-mounted assembly aid. The angled Multi-Unit abutments with rotation lock can be placed in four different positions, taking into account the orientation of the angulation. The straight Multi-Unit abutments without rotation lock can be placed freely in the implant. Once the correct seating of the Multi-Unit abutments in the implant has been ensured, the assembly aid can be removed by turning it anticlockwise and disposed of.

The assembly aid is not suitable for transmitting torque. Therefore, it must not be used to tighten the Multi-Unit abutments on the implant.

The tightening of the Multi-Unit abutments is done with the following OKTAGON instruments:

- Straight: insertion aid (REF 77703)
- Angled: Torx screwdriver (REF 63978, 64551, 64553, 77900)

The recommended torque of 35 Ncm must be observed. It is recommended to check the specified torque again after about a week and to tighten it if necessary.

Contraindications

In principle, general medical contraindications as well as local, absolute and relative contraindications for dental surgery must be observed.

Absolute contraindications include, among other things, serious internal illnesses, bleeding disorders, metabolic disorders, bone metabolism disorders (e.g. osteoporosis), reduced immune defence, pregnancy (and breastfeeding) as well as incomplete growth of the upper and lower jaw and blood clotting disorders. Furthermore, all components are made of materials that are suitable for the indicated applications. The use of the affected products is contraindicated in case of allergies or hypersensitivity to the chemical components of the following materials used:

- Titanium Grade 4 (3.7065)
- Titanium Grade 5 ELI (3.7165 Medical Grade, Ti-6Al-4V ELI, Grade 23)
- Titanium Grade 5 ELI (3.7165 Medical Grade, Ti-6Al-4V ELI, Grade 23) + Zirconium nitride (ZrN)

Relative contraindications are bruxism, psychiatric disorders, medication, drug, alcohol and tobacco use, as well as insufficient bone supply and poor bone quality, poor oral hygiene and untreated remaining teeth. In case of insufficient bone volume, we refer you to the common augmentation procedures and the corresponding literature. Incorrect selection and/or positioning of the prosthetic components can lead to biomechanical overloading due to unfavourable force distribution, resulting in implant loss, fatigue fractures of implants, abutments or prosthetic screws.

Notes

The user must ensure that all components used are secured against aspiration or swallowing.

All prosthetic components are designed for single-use. Reusing single-use products carries an increased risk of infection and risk-free functional safety (e.g. in terms of stability or accuracy of fit) cannot be guaranteed.

At every step of the prosthetic restoration (e.g. insertion of the cover screw or gingiva former, impression taking and insertion of the final prosthesis), it must be ensured that the internal geometry of the implant is free of contamination, so that the implant-secondary part connection and thus the success of the restoration are not negatively affected. All prosthetic components and instruments are supplied non-sterile and are labelled with the corresponding symbols. They must be cleaned, disinfected and sterilised before use. Only system-compatible components may be used. The colour coding and/or labelling must be observed. Furthermore, any kind of modification in the procedure or to the product can impair the success of the treatment or damage the product. The patient must be explicitly and thoroughly informed about the necessary oral hygiene measures. It is recommended that the patient be enrolled in a regular recall system. The OKTAGON prosthetic components have not been evaluated for safety and compatibility in the MRI environment and have not been tested for heating, surface diffusion or image distortion in the MRI environment. Safety is uncertain in the MRI environment. Scanning a patient with these products could result in patient injury.

Liability and warranty

The use of non-system components, as well as any kind of modification, can impair the function of the implant system (problems with fit, reduction in fatigue strength, as well as damage and even fracture of the implant) and excludes any guarantee or compensation by the manufacturer. This applies in particular to other application procedures that are not recommended. The processing and application of the products is beyond our control and is the sole responsibility of the user. Any liability for damage caused in this way is excluded. The technical advice on the use of our products is provided orally, in writing, by means of electronic media or by means of demonstrations. This is based on the state of the art in science and technology as known to us at the time the product is placed on the market. It does not release the user from the obligation to personally test the product for its indication-related suitability and to further educate themselves in the field of dentoalveolar surgery/implantology.

Any serious incident that has occurred in relation to the devices must be reported to the manufacturer and to the competent authority of the Member State in which the user and/or patient is established.

Disposal

When disposing of prosthetic components and instruments (at the end of their service life or after the specified shelf life has expired), care should be taken to ensure that the product is disposed of in the waste for biological hazardous substances. All packaging components are disposed of in accordance with national requirements (e.g. dual waste system).


Summary of safety and clinical performance

The 'Summary of Safety and Clinical Performance' (SSCP) provides a summary of the clinical data and other information on the safety and clinical performance of the OKTAGON prosthetic components. The SSCP can be viewed on the European Database on Medical Devices (EUDAMED) website <https://ec.europa.eu/tools/eudamed> after publication or otherwise requested by emailing RA@meisinger.de. The user is encouraged to inform the patient about the availability of the SSCP on the aforementioned website.

General information

Please also note the general application and safety instructions for MEISINGER products in the medical field and the instructions for reprocessing (cleaning, disinfection and sterilisation) of Hager & Meisinger GmbH medical devices.

Hager & Meisinger GmbH
Hansestr. 10
41468 Neuss • Germany
Tel. +49 2131 2012-0
Fax: +49 2131 2012-222
Email: info@meisinger.de
Internet: www.meisinger.de

Abbreviation	Significance
TL RP	Tissue Level Regular Platform
TL WP	Tissue Level Wide Platform
BL NC	Bone Level Narrow Connect
BL RC	Bone Level Regular Connect

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