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
Implant insertion

Following the final preparation of the implant site with the system's own instruments, the implantation of the selected implant can be performed in accordance with the manufacturer's information.

Determination of the drilling depth / selection of the stop sleeve

The required drilling depth and selection of the appropriate stop sleeve is determined before the start of the operation according to the selected implant length. With a required drilling depth of 11.0 mm for example, the stop sleeve CL035 must be selected. This has a length of 8.0 mm and allows a maximum drilling depth of 11.0 mm when mounted. The details of all maximum possible drilling depths of the stop sleeves included in the system can be found in the article overview under "Contents". The drilling depths are marked on the respective drill sleeve for clear identification. The length marking on the various drill bits is an additional optional orientation aid and should not primarily be used to determine the drilling depth. Important: The drilling depth marking on the 187RF 204 018 pilot core bit and the TDS0G 204 020 twist drill only starts at a depth of 8 mm.

ATTENTION: To avoid injuring anatomical structures, the apical overlength of the twist drills of up to 0.8 mm must be observed.

Preparation of the implant site

For the subsequent preparation of the implant site, the selected drill sleeve must be placed on the drill used in each step. Make sure that the drill sleeve is securely seated. In order to avoid overheating of the jawbone, sufficient cooling with preferably physiological, sterile saline solution must be ensured during preparation. In addition, the drills should be used intermittently and under low pressure.

Pilot drilling of the alveolar ridge

In order to prevent the twist drills from slipping when, for example, the alveolar ridge tapers to a point, the pilot drilling can be used beforehand.

Use of the extension drills

For the subsequent extension drilling, the twist drills are used in ascending diameter until the desired diameter of the implant site is reached. The first step is to drill a pilot bore to a diameter of 2.0 mm. This determines the alignment of the implant site in addition to the predefined depth.

Use of twist drills with guiding tip

Between the actual extension holes, intermediate holes are drilled using the twist drills with guiding tip. These drills have a non-toothed guiding tip, which enables a particularly safe extension of the implant site. When using the drills, initially only the crestal part of the implant site is extended, while the guiding tip uses the apical part of the drill hole as a guide channel. Once the desired diameter of the implant site is reached, the extension is completed with one of the fully toothed extension drills. The final diameter as well as the fine preparation of the implant site must be created with the instruments of the planned implant system. Depending on the existing case situation and the desired implant diameter, the instruments to be used can be selected individually. An exemplary sequence of instrument use is listed below (the stop sleeve to be used is not included):

- | | | |
|---|--|--|
| 1. Initial drilling:
• 187RF 204 018 | 3. Extension drilling, Ø 2.8 mm:
• TDS08 204 028
• TDS0G 204 028 | 5. Extension drilling, Ø 4.2 mm:
• TDS08 204 042
• TDS0G 204 042 |
| 2. Pilot drilling, Ø 2.0 mm:
• TDS0G 204 020 | 4. Extension drilling, Ø 3.5 mm:
• TDS08 204 035
• TDS0G 204 035 | |

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