Application and safety instructions for MEISINGER

Micro Screw T5 System

Micro Screw T5 System Basic

CONTENT

Language		Page
English (EN)		2-4
French (FR)		5-7
German (DE)		8-10
Italian (IT)		11-13
Spanish (ES)		14-16

Application and safety instructions for MEISINGER

Micro Screw T5 System Micro Screw T5 System Basic

Screw System for safe Stabilization of Cortical Bone Graft Transplants 'Olsberger Konzept' by Prof. Dr. Fouad Khoury

Art.-No.: MST10 4 – MST10 14, MST12 4 – MST12 14, BMST0 and BMSTB

The Bone Management Micro Screw T5 Systems are used for the safe fixation and stabilisation of cortical bone grafts. The special feature of these systems are the osteosynthesis screws, which have a diameter of 1.0 mm and 1.2 mm. These Micro Screws are made of surgical stainless steel, which gives them high stability despite their reduced diameter.

INDICATIONS FOR USE

The Bone Management Micro Screw T5 Systems are indicated for the treatment of bone fractures, especially for the fixation of transplanted bone blocks during the augmentation process in the oral cavity and maxillomandibular surgical field. Note: Micro Screws are not intended to remain in the body permanently. After they have fulfilled their supportive function such as is the case after healing of a transplant, or healing of a fracture, for example, they need to be removed completely.

APPLICATION

The recommended speeds for the application and the maximum speeds are summarized in the instrument overview under „Content“. For proper, safe use, these must be observed. The user must ensure that the screws are secured against aspiration both during the attachment and during its removal (e.g. by placing a rubber dam).

Placement of the pilot hole through bone cylinder and recipient bone

For the application of the Micro Screws, a pilot hole through bone graft and recipient bone is necessary.

This is done using a twist drill with a 0.2 mm smaller diameter:

- 1.0 mm Micro Screw $\pm$ twist drill MSPB1 or initial bur 202RF 008 ($\varnothing$ 0.8 mm)
- 1.2 mm Micro Screw $\pm$ twist drill MSPB2 or initial bur 202RF 010 ($\varnothing$ 1.0 mm)

The respective drilling should take place intermittently and under constant cooling with sterile, physiological saline solution. Since the Micro Screws have a self-tapping thread due to specially shaped thread flanks, no expansion drilling (gliding hole) to the outer diameter of the screws is necessary.

Selection of the screwdriver and picking up of the Micro Screws

Depending on the selected procedure (machine or manual screw insertion), the appropriate screwdriver must be selected. The screwdriver MSTDM is used for manual insertion. For machine insertion with the angle piece, the MSTDW is used.

Fixation of the bone graft

Using the selected screwdriver, a Micro Screw of appropriate length is picked up and the bone graft is attached to the recipient bone. The Micro Screw can be inserted into the bone without thread cutting, as it has a self-tapping thread design. The Micro Screw is tightened in a controlled manner until the bone cylinder is firmly fixed in the site. The bone must not be over stressed. In addition, during the insertion process, the axis specified by the pilot hole must be observed to prevent the screw from tilting. If a safe fixation of the graft with one screw is not possible, fixation with at least 2 screws is recommended to prevent any rotation.

By tilting the screwdriver sideways, the screwdriver can be released from the inserted Micro Screw. Do not pull it out vertically, as the strong bond between the Micro Screw and the screwdriver can lead to the Micro Screw being pulled out.

Removing the Micro Screws

After clinically confirmed healing and functional restoration, the Micro Screws must be removed using a screwdriver from the Micro Screw T5 System in an anticlockwise direction. If screw loosening, infection or weakening of the bone tissue in the area of the inserted bone screws occurs, the Micro Screws must be removed. The decision on the time of removal is the responsibility of the treating physician.

Sterilisation of the Micro Screws

The instruments and Micro Screws are supplied non-sterile and must be prepared (cleaning, disinfection and sterilization) prior to clinical use. For proper processing of the items, please follow the instructions for processing (cleaning, disinfection and sterilization) of medical devices from Hager & Meisinger GmbH at <https://www.meisinger.de/services/downloads>.

Disposal

Ensure that the product is disposed of with biohazardous waste at the end of use or at the end of the specified shelf life. All components of the packaging must be disposed of in accordance with national regulations (e.g. the dual waste system).

Intended Use

The products for bone fixation are intended for the fixation and stabilization of bone segments (all screws) and for the support of barrier membranes (Tenting-Screws) in dental surgery. The screws are used for the restoration of anatomical structures and for pre-implant vertical or horizontal augmentation.

Intended Patient Population

Taking contraindications into account, the products can be used for patients of any gender, when the alveolar arch is developmentally mature.

Intended User Profile

Professional users / specialist users in the field of dental implantology.

Contraindications

Basically, general medical as well as local, absolute and relative contraindications for dental surgical procedures must be considered.

Absolute contraindications

- Active infections in or near the area to be augmented as well as local or systemic pathological processes (e.g. symptoms such as fever, local inflammation, abscesses)
- Insufficient bone volume (quality / quantity)
- Diseases or medication that impair bone metabolism

Relative contraindications

- Dentoalveolar growth which has not come to an end (exception: cases in which no dentoalveolar growth can be expected, e.g. ectodermal dysplasia)
- Known allergies and/or foreign body reaction to the alloy composition of the Micro Screws, such as nickel
- Systemic and/or metabolic diseases or medical treatments that lead to progressive deterioration of the bone (e.g. cortisone, immunosuppressants, bisphosphonates, diabetes, etc.)
- Drug and alcohol abuse
- Lack of patient compliance
- Poor circulation
- Hard physical labor or active sports
- Mental condition, which can lead to disregard of the medical order
- Highly atrophic jaw

Attention must be paid to the protection of the anatomical structures (safety distance min. 2 mm) as well as to the course of the adjacent teeth / tooth roots (risk of damage, infection, resorption). To avoid infection and resorption, care should be taken to avoid contact between the graft and adjacent teeth. Mechanical stability is essential for successful bone healing..

SEPARATE INSTRUCTIONS FOR SCREWS: The bone screws are compatible with the instruments of the Micro Screw T5 System and Micro Screw T5 System Basic. The screws may only be used with these system instruments. When selecting the screw length as well as the screw placement, the evaluation of the available bone by means of X-ray diagnostics must be taken into account. The screw head should be countersunk as flat as possible in the bone block to prevent soft tissue perforation or dehiscence. Excessive physical activity as well as trauma impairing the augmentation site can lead to premature failure of the screws due to loosening or breaking. The screws may only be used in accordance with the application and safety instructions. Any kind of modification in the procedure or on the product can lead to impairment of the success of the treatment or damage to the product. The screws are a disposable product and must therefore not be reprocessed and reused after use. Even with screws that have already been used but are visually undamaged, the previous mechanical stress can cause microscopic damage, which can lead to fatigue and product failure. The Micro Screws have not been evaluated for safety and compatibility in the MR environment. They have not been tested for heating, migration, or image artifact in the MR environment. The safety of the Micro Screws in the MR environment is unknown. Scanning a patient who has this medical device may result in patient injury.

The user must ensure that the pins are secured against aspiration both during the attachment of the membrane and during the removal of the pins.

The 'Summary of Safety and Clinical Performance' (SSCP) summarizes the clinical data and other information on the safety and clinical performance of the pins. On the website of the European Database on Medical Devices (EUDAMED) <https://ec.europa.eu/tools/eudamed> the SSCP can be viewed after publication or otherwise requested by sending an e-mail to RA@meisinger.de. The user is obliged to inform the patient of the availability of the SSCP on the aforementioned website.

Serious incidents related to the medical device shall be notified to the manufacturer and to the competent authority of the Member State where the user and / or patient is located.

General instructions

Please follow general application and safety instructions for medical devices produced by MEISINGER and also the advice for processing (cleaning, disinfection and sterilization) of medical devices from Hager & Meisinger GmbH at <https://www.meisinger.de/services/downloads>.

REF

MD

LOT

Do not reuse

Non-sterile

Caution

CE / CE 0044

Do not use if the packaging is damaged

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
* Minimal diameter External diameter

** Surgical stainless steel according to DIN ISO 5832-1, ASTM F138. For further information on alloy composition please contact the manufacturer.

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