Application and safety instructions for MEISINGER
Micro Screw T5 System Plus

CONTENT

Language		Page
English (EN)		2-4
French (FR)		5-7
German (DE)		8-10
Italian (IT)		11-13
Spanish (ES)		14-16

Application and safety instructions for MEISINGER

Micro Screw T5 System Plus

Screw System for safe Stabilization of Cortical Bone Graft Transplants 'Olsberger Konzept' by Prof. Dr. Fouad Khoury

Art.-No.: BMTSP

The Micro Screw T5 System Plus is a combined system of the Trephine Ejection Kit, which enables the minimally invasive extraction of local bone cylinders, and the Micro Screw T5 System, which is used for the safe fixation and stabilization of bone grafts. The Micro Screw T5 System Plus thus enables the complete workflow of the Bone-Core-Technique to be carried out. For the extraction of local bone cylinders, the system contains two-part, internally cooled trephines, which, thanks to the removable trephine working parts, enable the bone cylinders to be pushed out particularly easily with the help of special ejection instruments. For the subsequent fixation of the bone cylinders, the Micro Screw T5 System Plus offers osteosynthesis screws with a diameter of 1.0 mm or 1.2 mm as well as all the necessary instruments for setting inserting the screws. The osteosynthesis screws are made of surgical stainless steel, which gives them high stability despite their reduced diameter. The combination of the two techniques in one system and the optimally matched instruments offer the user a holistic, clear workflow when performing the Bone-Core-Technique and allow for the highest precision from the extraction to the fixation of autologous bone cylinders.

INDICATIONS FOR USE

The Bone Management Micro Screw T5 System Plus is indicated for the intraoral extraction of bone cylinders and the treatment of bone fractures, especially for the fixation of transplanted bone blocks during the augmentation process in the oral cavity and maxillomandibular surgical field. Note: Micro Screws are not intended to remain in the body permanently. After they have fulfilled their supportive function e. g. after a transplant or fracture has healed, they need to be removed completely.

APPLICATION

The recommended speeds for the application and the maximum speeds are summarized in the instrument overview under „Content“. For proper and safe use, these must be observed. The user must ensure that the screws are secured against aspiration both during the attachment and during its removal (e. g. by placing a rubber dam).

Before application


The instruments and Micro Screws are non-sterile when delivered and must be processed (cleaned, disinfected and sterilized) according to validated instructions for processing before use.

For cleaning and disinfection, the products are disassembled as far as possible. For that the shank of the trephine is attached to the hand wheel, so the trephine working part can be released and removed by a left-hand turn. The secure fit of the working part has to be checked. Subsequently, the instruments can be sterilized according to the instructions for processing. The secure fit of the working part has to be checked.

Since the trephines have an internal cooling, these are instruments with increased processing requirements. According to the RKI guideline, mechanical cleaning and disinfection in disinfection devices and sterilization with moist heat is recommended. For proper processing of the articles, please note the instructions for processing (cleaning, disinfection and sterilization) of medical devices produced by Hager & Meisinger GmbH.

Initial Drilling

According to the planned diameter of the bone cylinder the appropriate initial drill is selected. With this drill the bone at the donor site is punch-marked in order to avoid slippage of the trephine.


Extraction of the bone cylinder from the donor site

According to the planned diameter of the bone cylinder the appropriate trephine is selected. Before use, the fixed position of the trephine working part on the shank has to be checked.

Caution: Depending on the construction, the trephines may not be used counterclockwise.


At the punch-marked donor site the bone cylinder of the planned length can now be extracted with the aid of the trephine. For this purpose, the trephine is placed on the cortical bone and intermittently inserted into the bone. It is recommended to support the working hand to control the trephine optimally. The speed should not be too low to avoid tilting of the trephine in the depth of the bone. Please note the recommended and maximum speeds. The depth marks on the trephine working parts serve for length orientation.

Extraction of the bone cylinder from the trephine

To extract the bone cylinder from the trephine working part, the trephine is removed from the handpiece. To achieve a better handling, the ejection needle is attached to the hand wheel. Then, it can be led through the lumen of the trephines' internal cooling from the shank end in the direction of the trephine working part, so the bone cylinder is pushed out.

Caution: In the case of soft bone, the application of the ejection needle is contraindicated due to its small diameter. In this case, proceed according to step "Alternative Extraction of the bone cylinder from the trephine in case of soft bone".

Alternative Extraction of the bone cylinder from the trephine in case of soft bone

In case of soft bone, the extraction of the bone cylinder is performed with the aid of the ejection instruments. For that, the shank of the trephine is attached to the hand wheel, so that the trephine working part can be released and removed by a left-hand turn. Then, the trephine working part is inserted into the ejection sleeve. With the aid of the appropriate ejection instrument the bone cylinder is pushed out.

The 2.1 mm and 2.5 mm diameter trephines go with the 2.1 mm diameter ejection instrument.

The 2.9 mm and 3.1 mm diameter trephines go with the 2.9 mm diameter ejection instrument.

Placement of the pilot hole through bone cylinder and recipient bone

For the application of the Micro Screws, a pilot hole through bone graft and recipient bone is necessary.

This is done using a twist drill with a 0.2 mm smaller diameter:

- 1.0 mm Micro Screw $\triangleq$ twist drill MSPB1 or initial bur 202RF 008 ($\varnothing$ 0.8 mm)
- 1.2 mm Micro Screw $\triangleq$ twist drill MSPB2 or initial bur 202RF 010 ($\varnothing$ 1.0 mm)

The respective drilling should take place intermittently and under constant cooling with sterile, physiological saline solution. Since the Micro Screws have a self-tapping thread due to specially shaped thread flanks, no expansion drilling (gliding hole) to the outer diameter of the screws is necessary.

Selection of the screwdriver

Depending on the selected procedure (machine or manual screw insertion), the appropriate screwdriver must be selected. The screwdriver MSTDM is used for manual insertion. For machine insertion with the angle piece, the MSTDW is used.

Fixation of the bone graft

Using the selected screwdriver, a Micro Screw of appropriate length is picked up and the bone graft is attached to the recipient bone. The Micro Screw can be inserted into the bone without thread cutting, as it has a self-tapping thread design. The Micro Screw is tightened in a controlled manner until the bone cylinder is firmly fixed in the site. The bone must not be over stressed. In addition, during the insertion process, the axis specified by the pilot hole must be observed to prevent the screw from tilting. If a safe fixation of the graft with one screw is not possible, fixation with at least 2 screws is recommended to prevent any rotation.

By tilting the screwdriver sideways, the screwdriver can be released from the inserted Micro Screw. Do not pull it out vertically, as the strong bond between the Micro Screw and the screwdriver can lead to the Micro Screw being pulled out.

After use

The instruments need to be cleaned, disinfected and sterilized after every application. For cleaning and disinfection, the instruments are disassembled as far as possible. With the aid of a tweezer the rubber rings on the shank are now removed and disposed. Then, the instruments can be reprocessed.

After completed cleaning and disinfection, preferably completely mechanical, a new rubber ring (Art.-No. 2151) is mounted on the shank and the trephine working part is again attached and fixed to the shank by a right-hand turn. Then, the instruments can be sterilized according to the instructions for processing. For a proper processing of the articles please note the instructions for processing (cleaning, disinfection and sterilization) of medical devices from Hager & Meisinger GmbH.

Removing the Micro Screws

After clinically confirmed healing and functional restoration, the Micro Screws must be removed using a screwdriver from the Micro Screw T5 System in an anticlockwise direction. If screw loosening, infection or weakening of the bone tissue in the area of the inserted bone screws occurs, the Micro Screws must be removed. The decision on the time of removal is the responsibility of the treating physician.

Disposal

Ensure that the product is disposed of with biohazardous waste at the end of use or at the end of the specified shelf life. All components of the packaging must be disposed of in accordance with national regulations (e.g. the dual waste system).

Intended Use

The products are intended for the intraoral extraction of bone cylinders and the fixation and stabilization of bone segments in dental surgery. The screws are used for the restoration of anatomical structures and for pre-implant vertical or horizontal augmentation.

Intended Patient Population

Taking contraindications into account, the products can be used for patients of any gender, when the alveolar arch is developmentally mature.

Intended User Profile

Professional users / specialist users in the field of dental implantology.

Contraindications

Basically, general medical as well as local, absolute and relative contraindications for dental surgical procedures must be considered.

Absolute contraindications

- Active infections in or near the area to be augmented as well as local or systemic pathological processes (e.g. symptoms such as fever, local inflammation, abscesses)
- Insufficient bone volume (quality / quantity)
- Diseases or medication that impair bone metabolism

Relative contraindications

- Dentoalveolar growth which has not come to an end (exception: cases in which no dentoalveolar growth can be expected, e.g. ectodermal dysplasia)
- Known allergies and/or foreign body reaction to the alloy composition of the Micro Screws, such as nickel
- Systemic and/or metabolic diseases or medical treatments that lead to progressive deterioration of the bone (e.g. cortisone, immunosuppressants, bisphosphonates, diabetes, etc.)
- Drug and alcohol abuse
- Lack of patient compliance
- Poor circulation
- Hard physical labor or active sports
- Mental condition, which can lead to disregard of the medical order
- Highly atrophic jaw

Attention must be paid to the protection of the anatomical structures (safety distance min. 2 mm) as well as to the course of the adjacent teeth / tooth roots (risk of damage, infection, resorption). To avoid infection and resorption, care should be taken to avoid contact between the graft and adjacent teeth. Mechanical stability is essential for successful bone healing. For all cutting and drilling procedures (transplant removal / other work), the following applies: In order to reduce the risk of overheating the bone and necrosis formation along with it, the respective instrument has to be used intermittently using little pressure and continuous cooling using a sterile physiological saline solution. It should be avoided to trepan the lingual cortical bone as damage to the soft tissue and vessels of the mouth base can lead to serious complications. The extracted bone graft should be used as quickly as possible in the recipient site or stored in blood or isotonic saline solution. It is recommended to fill larger defects after bone extraction. In the donor site, a tension-free and complete wound closure must follow quickly.

SEPARATE INSTRUCTIONS FOR SCREWS: The bone screws are compatible with the instruments of the Micro Screw T5 System Plus. The screws may only be used with these system instruments. When selecting the screw length as well as the screw placement, the evaluation of the available bone by means of X-ray diagnostics must be taken into account. The screw head should be countersunk as flat as possible in the bone block to prevent soft tissue perforation or dehiscence. Excessive physical activity as well as trauma impairing the augmentation site can lead to premature failure of the screws due to loosening or breaking. The screws may only be used in accordance with the application and safety instructions. Any kind of modification in the procedure or on the product can lead to impairment of the success of the treatment or damage to the product. The screws are a disposable product and must therefore not be reprocessed and reused after use. Even with screws that have already been used but are visually undamaged, the previous mechanical stress can cause microscopic damage, which can lead to fatigue and product failure. The Micro Screws have not been evaluated for safety and compatibility in the MR environment. They have not been tested for heating, migration, or image artifact in the MR environment. The safety of the Micro Screws in the MR environment is unknown. Scanning a patient who has this medical device may result in patient injury.

The user must ensure that the screws are secured against aspiration both during the attachment of the membrane and during the removal of the screws.

The 'Summary of Safety and Clinical Performance' (SSCP) summarizes the clinical data and other information on the safety and clinical performance of the products. On the website of the European Database on Medical Devices (EUDAMED) <https://ec.europa.eu/tools/eudamed> the SSCP can be viewed after publication or otherwise requested by sending an e-mail to RA@meisinger.de. The user is obliged to inform the patient of the availability of the SSCP on the aforementioned website.

Serious incidents related to the medical device shall be notified to the manufacturer and to the competent authority of the Member State where the user and / or patient is located.

General instructions

Please follow general application and safety instructions for medical devices produced by MEISINGER and also the advice for processing (cleaning, disinfection and sterilization) of medical devices from Hager & Meisinger GmbH at <https://www.meisinger.de/services/downloads>.

REF

MD

LOT

Do not reuse

Non-sterile

Caution

CE / CE 0044

Conformity marks

Do not use if the packaging is damaged

Date of manufacture

Packaging unit

Consult instructions for use

ifu.meisinger.de

Hager & Meisinger GmbH

Hansemannstr. 10

41468 Neuss • Germany

Tel. +49 2131 2012-0

Fax: +49 2131 2012-222

Email: info@meisinger.de

Internet: www.meisinger.de

UDI

Unique Device Identification

UK RP

Swiss authorised representative

CH REP

UK authorised representative

Content

	Initial Bur		Twist Drill		Screwdriver RA	Screwdriver manual
						
	1	1	1	1	1	1
Fig.	202RF	202RF	MSPB1	MSPB2	MSTDW	MSTDW
Shank ¹	206	206	204	204	204	-
Size	008	010	008	010	-	-
Length mm	14.0	14.0	14.0	14.0	18.0	148.0
	-	-	-	-	-	-
	0.8	1.0	0.8	1.0	-	-
Opt. speed rpm	1000	1000	1000	1000	10	-
Max. speed rpm	6000	6000	6000	6000	15	-
MD	CE 0044					CE

	Screws* **					
						
	4	4	4	4	4	4
Fig.	MST10 4	MST10 6	MST10 8	MST10 10	MST10 12	MST10 14
Length mm	4.0	6.0	8.0	10.0	12.0	14.0
	0.6	0.6	0.6	0.6	0.6	0.6
	1.0	1.0	1.0	1.0	1.0	1.0
Head-Ø mm	1.7	1.7	1.7	1.7	1.7	1.7
MD	CE 0044					

	Screws* **					
						
	1	1	1	1	1	1
Fig.	MST12 4	MST12 6	MST12 8	MST12 10	MST12 12	MST12 14
Length mm	4.0	6.0	8.0	10.0	12.0	14.0
	0.7	0.7	0.7	0.7	0.7	0.7
	1.2	1.2	1.2	1.2	1.2	1.2
Head-Ø mm	1.7	1.7	1.7	1.7	1.7	1.7
MD	CE 0044					

¹ 204 = RA

Minimal diameter




External diameter

** Surgical stainless steel according to DIN ISO 5832-1, ASTM F138. For further information on alloy composition please contact the manufacturer.

Content

	Initial bur trephines				Shank***
					
	1	1	1	1	4
Fig.	230KH	230KH	230KH	230KH	229KH
Shank ¹	204	204	204	204	204
Size ²	021	025	029	033	000
Length mm	2.0	2.0	2.0	2.0	20.5
	2.1	2.5	2.9	3.3	-
	3.1	3.5	3.9	4.3	-
Opt. speed rpm	300-500	300-500	300-500	300-500	300-500
Max. speed rpm	800	800	800	800	4000
MD	CE 0044				

	Trephines			
				
	1	1	1	1
Fig.	229KH	229KH	229KH	290KH
Shank ¹	900	900	900	900
Size ²	021	025	029	033
Length mm	15.5	15.5	15.5	15.5
	2.1	2.5	2.9	3.3
	3.1	3.5	3.9	4.3
Opt. speed rpm	300-500	300-500	300-500	300-500
Max. speed rpm	2000	2000	2000	2000
MD	CE 0044			

	Ejection instruments			Ejection sleeve	Hand wheel
					
	1	1	1	1	1
Fig.	ASTKH	ASTKH	ASNKH	ASHKH	HRKH1
Size ²	021	029	008	-	-
Length mm	22.0	22.0	37.0	28.0	15.5
	-	-	-	5.5	-
	5.5	5.5	-	-	-
	2.1	2.9	0.8	-	-
MD	CE				

¹ 204 = RA ² Largest working part diameter in 1/10 mm

 Minimal diameter  Internal diameter  External diameter

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