Application and safety instructions for MEISINGER

Tenting Screw System Umbrella Screw

85FLBM39 – Rev. 2 10/25

CONTENT

Language		Page
English (EN)		2-4
French (FR)		5-7
German (DE)		8-10
Italian (IT)		11-13
Spanish (ES)		14-16

Application and safety instructions for MEISINGER

Tenting Screw System Umbrella Screw

Bone Screws

Art.-No.: BTS00 and USL600

The Bone Management Tenting Screw System and Umbrella Screw are a treatment options for both simple and complex bone defects and offers an alternative to the shell technique and the use of titanium-reinforced membranes. The screws contained in the System as well as the Umbrella Screw are made of the titanium alloy Grade 5-ELI (Titanium Grade 23) and have a self-tapping thread. The screw head includes an interface to accommodate the insertion tool (hexalobular socket).

INDICATIONS FOR USE

The bone screws are indicated for the fixation and stabilization of bone grafts and bone substitutes as well as for the support of barrier membranes. They are used in dental surgery therapies (osteosynthesis, correction of degenerative diseases) to restore anatomical structures, but are not suitable as a replacement for such structures.

UTILIZATION

The bone screws can be used for a variety of surgical techniques depending on the case situation. The selection of the technique and the composition of the instruments to be used must be selected and prepared for use before the operation begins. The user must ensure that the screws are secured against Aspiration both during fixation and removal (e.g. by placing a rubber dam). With the Umbrella Screw, it is also, for example, possible to secure the screw against aspiration using suture material.

Preparation of the screws

Due to the self-tapping thread design, no pilot hole is usually required for inserting the screws. If required, the bone can be pre-punched using the initial bur included in the system. Pre-graining is recommended for hard bone. In the case of very hard bone, pre-drilling using the drill included in the system may be indicated.

Selecting the screwdriver and holding the screws

Depending on the selected procedure (mechanical or manual screw insertion), the appropriate screwdriver from the manufacturer must be selected. The TL0T1 screwdriver is used for manual insertion. The 31053 screwdriver is used for automatic insertion with the contra-angle handpiece. The use of third-party accessories may impair the function of the screws and invalidates any warranty or replacement.

The screws in the Tenting Screw System are contained in a sterilizable box that can be opened using a sliding mechanism. The screws are removed from the box using the screwdriver. To do this, the screwdriver is inserted vertically into the interface contained in the screw heads. The screws are not fixed in the box. Therefore, care must be taken to ensure vibration-free working both when opening the box and when picking up the screws with the screwdriver.

As the screws are implantable products, the screws must only be picked up using the screwdriver to avoid touching them with your hand.

Postoperative procedure

It is recommended to remove the screws counterclockwise using a suitable screwdriver after clinically confirmed healing and functional restoration. If an infection or weakening of the bone tissue occurs in the area of the inserted bone screws, the screws must be removed. Explantation should also be carried out in the event of screw loosening and generally in children to avoid growth inhibition. The decision on removal and the time of removal is the responsibility of the attending physician

Sterilisation of the screws

The screws are supplied non-sterile and must be processed before use (cleaning / disinfection / sterilization). For proper processing of the items, please observe the instructions for processing (cleaning, disinfection and sterilization) of medical devices from Hager & Meisinger GmbH.

Disposal

When disposing of the instruments (at the end of their service life or after expiry of the specified shelf life), ensure that the product is disposed of in the waste for biohazardous substances. All packaging components are disposed of in accordance with national regulations (e.g. dual waste system).

Intended Use

The bone fixation products are intended for the fixation and stabilization of bone segments (all screws) and for supporting membranes (tenting and umbrella screws) in dental surgery. The screws are used to restore anatomical structures and for vertical or horizontal pre-implant augmentation.

Intended Patient Population

Taking contraindications into account, the products can be used for patients of any gender, when the alveolar arch is developmentally mature.

Intended User Profile

Professional users / specialist users in the field of dental implantology.

Contraindications

In principle, general medical as well as local, absolute and relative contraindications to dental surgery must be considered.

Absolute contraindications

- Active infections in or near the area to be augmented as well as local or systemic pathological processes (e.g., symptoms such as fever, local inflammation, abscesses)
- Insufficient bone supply (quality / quantity)
- Illnesses or medication that impair bone metabolism
- Insufficient soft tissue supply, which does not allow complete coverage of the augmented area

Relative contraindications

- Dentoalveolar growth not completed (exception: cases where no dentoalveolar growth is expected, e.g., ectodermal dysplasia) cases where no dentoalveolar growth is expected, e.g., ectodermal dysplasia)
- Known allergies and / or foreign body reaction to the implant material
- Systemic and/or metabolic illnesses or medical treatments that lead to progressive bone deterioration (e.g., cortisone, immunosuppressants, bisphosphonates and diabetes)
- Drug and alcohol abuse
- Lack of patient compliance
- Poor circulation
- Heavy work or active Sports
- Mental state which may lead to non-compliance with the medical order
- Highly atrophied jaw
- Nicotine consumption, which may jeopardize the success of the procedure by delaying bone wound healing

IMPORTANT: Care must be taken to protect the anatomical structures (safety distance min. 2 mm) as well as the course of the adjacent teeth / tooth roots (hazard of damage, infection / resorption). To avoid the development of infections and dehiscence, care should be taken to prevent contact between the graft and adjacent teeth. Cases in which the mechanical stability and mechanical rest of the graft cannot be achieved may compromise successful bone healing.

SEPERATE INSTRUCTIONS FOR SCREWS: When selecting the screw length as well as the screw placement, the evaluation of the available bone supply by means of X-ray diagnostics must be taken into account. To prophylactically counteract soft tissue perforation or dehiscence, tension-free and complete wound closure must be achieved. The entire augmentation area, including the bone screws, must not be loaded during the healing process. Failure to do so may result in bending, loosening, loss and/or breakage of the screws and this means a negative impact on the augmentation result. Also, excessive physical activity as well as trauma affecting the augmentation site can lead to premature failure of the screws by loosening or breaking. The screws must only be used in accordance with the instructions for use and safety. The complete sinking of the screw can be limited by very hard bone quality. Any kind of modification in the procedure or on the product can lead to impairment of the treatment success or damage to the product. The screws are a single-use product and must therefore not be reprocessed and reused after use. Even screws that have already been used but are visually undamaged may suffer microscopic damage from previous mechanical stress leading to fatigue and product failure. The screws have not been evaluated for safety and compatibility in the MR environment. They were not tested for heating, migration, or image artefacts in the MR environment. The safety of the screws is uncertain in the MR environment. Scanning a Patient with this medical device could result in injury to the patient.

Any serious-incident that has occurred in relation to the device should be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

The 'Summary of Safety and Clinical Performance' (SSCP) summarizes the clinical data and other information on the safety and clinical performance of the product. On the website of the European Database on Medical Devices (EUDAMED) <https://ec.europa.eu/tools/eudamed> the SSCP can be viewed after publication or otherwise requested by sending an e-mail to RA@meisinger.de. The user is obliged to inform the patient of the availability of the SSCP on the aforementioned website.

General instructions

Please follow general application and safety instructions for medical devices produced by MEISINGER and also the advice for processing (cleaning, disinfection and sterilization) of medical devices from Hager & Meisinger GmbH at <https://www.meisinger.de/services/downloads>.

Manufacturer

Catalogue number

Medical device

Batch Code


Do not reuse

Non-sterile

Caution

Conformity marks

Do not use if the packaging is damaged

Date of manufacture

Packaging unit

Consult instructions for use

Hager & Meisinger GmbH

Hansemannstr. 10

41468 Neuss • Germany

Tel. +49 2131 2012-0

Email: info@meisinger.de

Internet: www.meisinger.de

Unique Device Identification

Swiss authorized representative

Content BTS00

Screws


	7	7	7
Fig.	TS035 040 030	TS035 040 040	TS035 040 050
L ₁ (mm)	3.5	3.5	3.5
L ₂ (mm)	4.0	4.0	4.0
L ₃ (mm)	3.0	4.0	5.0
L ₄ (mm)	7.75	8.75	9.75
	CE 0044 /		

Screws


	3	3	3	2	2	2
Fig.	US050 040 070	US050 040 090	US050 040 110	US070 040 070	US070 040 090	US070 040 110
L ₁ (mm)	5.0	5.0	5.0	7.0	7.0	7.0
L ₂ (mm)	4.0	4.0	4.0	4.0	4.0	4.0
L ₃ (mm)	2.25	4.25	6.25	2.25	4.25	6.25
L ₄ (mm)	7.0	9.0	11.0	7.0	9.0	11.0
	CE 0044 /					

	1	1	1
Fig.	186RF	203RF	31053
Schank ¹	204	205	205
Size ²	018	011	-
Length mm	12.0	9.0	27.0
	-	-	-
	1.8	1.1	-
Opt. Speed rpm	800	1000	10
Max. Speed rpm	6000	6000	15
	CE 0044		

Screwdriver


	1
Fig.	TL0T1
Length mm	148,0
	CE

USL600

Umbrella Screw

	5
Fig.	USL600
Length mm	9.0
Head-Ø mm	6.0
	1.4
	CE 0044 /

¹ 204 = RA, 205 = RA L ² Largest working part diameter in 1/10 mm

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