Application and safety instructions for the MEISINGER Master-Pin-X-Control Systems

CONTENT

Language	Page
English (EN)	2-4
French (FR)	5-7
German (DE)	8-10
Italian (IT)	11-13
Spanish (ES)	14-16

BMPXB

BMPXU

Application and safety instructions for MEISINGER Master-Pin-X-Control Systems

Pin systems for membrane fixation

Art.-No.: BMPXB, BMPXU

Master-Pin-X-Control Systemes has been developed especially for the fixation of resorbable, non-resorbable and titanium-reinforced membranes. The pins are made of grade 5 titanium alloy and have a self-cutting thread that allows them to be inserted easily and precisely into the cortical bone.

APPLICATION INSTRUCTIONS

The Bone Management® Master-Pin-X-Control Systems are indicated for the fixation of the apical membrane ends to the local bone in order to avoid micromobility of the membrane. It enables the fixation of resorbable, non-resorbable and titanium-reinforced membranes.

APPLICATION

The speeds recommended for the application and the maximum permitted speeds are summarized in the article overview under "Contents". These must be observed to ensure proper and safe use.

Setting bleeding points

To freshen the bone, bleeding points can be placed in the cortical bone in the region to be augmented using the included twist drills.

Application of the augmentation material

The augmentation material is applied to the defect in accordance with the manufacturer's instructions and the selected membrane is adapted to the defect.

Recording the pins

The pins are contained in a sterilizable box that can be opened using a sliding mechanism. The pins are removed from the box using the TL0T1 screwdriver. In doing so, it is important to ensure a secure retention between the screwdriver and the pin. The pins are not fixed in the box. For this reason, care must be taken to ensure vibration-free work both when opening the box and when picking up the pins with the screwdriver.

⚠ As the pins are implantable products, they should be picked up only using the screwdriver in order to avoid touching them by hand.

Securing the membrane

To secure the membrane over the defect, the pins are applied through the membrane into the cortical bone using the TL0T1 screwdriver. In doing so, the pin can be inserted into the bone without pre-drilling in order to prevent aspiration, as it has a self-tapping thread design. The pin is gently tightened until the membrane is firmly fixed. The bone must not be overloaded in the process. To make it easier to remove the screwdriver from the pin after screwing it in, the screwdriver can be angled at 45°.

⚠ Compression of the membrane must be avoided.

Removing the pins

After a clinically confirmed, successful healing phase, a small incision is made in the gingiva above the pin to be removed. The pin can then be carefully loosened using the TL0T1 screwdriver by turning it counterclockwise. Depending on the case situation and used, the pins can be left in situ. If infection or weakening of the bone tissue occurs due to the inserted pins, the pins must be removed. Explantation is also essential if the pins become loose.

⚠ The pins are a disposable product and should therefore not be reused. The membrane must be handled in accordance with the manufacturer's instructions.

Sterilization of the pins

The instruments and pins are supplied non-sterile and must be processed before use (cleaning / disinfection / sterilization). For proper processing of the items, please observe the instructions for processing (cleaning, disinfection and sterilization) of medical devices from Hager & Meisinger GmbH.

Disposal

When disposing of the instruments (at the end of their service life or after expiry of the specified shelf life), ensure that the product is disposed of in the waste for biohazardous substances. All packaging components are disposed of in accordance with national regulations (e.g. dual waste system).

Contraindications

In principle, general medical as well as local, absolute and relative contraindications for dental surgery must be taken into account.

Absolute contraindications

- Incomplete dentoalveolar growth (exception: cases where no dentoalveolar growth is to be expected, e.g. ectodermal dysplasia)
- Acute infections in or near the area to be augmented as well as local pathological processes (e.g. symptoms such as fever, local inflammation, abscesses)
- Insufficient bone quality and/or quantity indications

Relative contraindications

- Known allergies and / or foreign body reaction to the implant material
- Systemic diseases
- Diseases that impair bone metabolism
- Medical treatments that can affect bone quality and quantity (e.g. cortisone, immunosuppressants)
- Drug and alcohol abuse
- Lack of cooperation from the patient
- Poor blood circulation
- Heavy work or active sports
- Mental condition that may lead to non-compliance with the doctor's orders
- Highly atrophied jaw

IMPORTANT: Care must be taken to protect the anatomical structures (safety distance min. 2 mm) and the course of the neighboring teeth / tooth roots (risk of damage, infection / dehiscence).

SPECIAL NOTES FOR PINS: Excessive physical activity and trauma that affects the membrane can lead to premature failure of the pin due to loosening or breakage. The pins may only be used in accordance with the application and safety instructions. Any type of modification to the procedure or the product may impair the success of the treatment or damage the product. The pins are a single-use product and after being used, must therefore not be reprocessed or reused. Even if the pins have already been used but are visually undamaged, the previous mechanical stress can lead to microscopic damage and thus to signs of fatigue and product failure. The pins were not evaluated for safety and compatibility in the MR environment. They have not been tested for heating, migration or image artifacts in the MR environment. The safety of the pins in the MR environment is uncertain. Scanning a patient with this medical device could result in injury to the patient.

RISK OF ASPIRATION: The user must ensure that the pins are secured against aspiration both when attaching the membrane and when removing it.

NOTE ON MEMBRANE TECHNOLOGY: Please refer to the information provided by the manufacturer of the membrane you have selected.

Serious incidents related to the device shall be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

The „Summary of Safety and Clinical Performance“ (SSCP) summarizes the clinical data and further information on the safety and clinical performance of the screws and can be provided as a copy upon request by e-mail to RA@meisinger.de. Once the European Database for Medical Devices (EUDAMED) is operational, the SSCP for the screws can be viewed on the following website <https://ec.europa.eu/tools/eudamed>.

General information

Please also observe the general application and safety instructions for MEISINGER products in the medical field and the instructions for processing (cleaning, disinfection and sterilization) of medical devices from Hager & Meisinger GmbH.

Manufacturer

Item number

Medical device

Batch number

Do not re-use

Non-sterile

Under U.S. federal law, this device may only be sold to or on behalf of trained medical professionals.

Caution

Conformity mark

Do not use if the packaging is damaged

Manufacturing date

Packaging unit

Follow the electronic instructions for use

ifu.meisinger.de

Hager & Meisinger GmbH
Hansemannstr. 10
41468 Neuss • Germany
Tel. +49 2131 2012-0
Fax: +49 2131 2012-222
Email: info@meisinger.de
Internet: www.meisinger.de


CONTENT

Master-Pin-X-Control Basic

Art.-No. BMPXB 
1x TL0T1 10x MP30
1x 203RF 006
1x 203RF 008
2x 203S 012

Master-Pin-X-Control Upgrade

Art.-No. BMPXU 
1x TL0T1 20x MP30

	Twist drill		Twist drill with stop	Pins with thread
				
Fig.	203RF	203RF	203S	MP30*
Shank ¹	205	205	204	-
Size ²	006	008	012	-
Length mm	7.0	7.0	4.0	3.20
	-	-	-	-
	0.6	0.8	1.2	1.45
Optimal speed rpm	1000	1000	1000	-
Max. Speed rpm	6000	6000	6000	-
CE 0044				

	Screwdriver
	
Fig.	TL0T1
	Figure 1:2 (length mm: 148.0)
	CE

*  ¹ 204=RA, 205=RA L ² Largest working part diameter in 1/10 mm  Minimum diameter  Outer diameter

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