Master-Pin-X-Control Plus

CONTENT

Language	Page
English (EN)	2 - 7
French (FR)	8 - 13
German (DE)	14 - 19
Italien (IT).....	20 - 25
Spanish (ES).....	26 - 31

English (EN)

Application and safety instructions for MEISINGER

Master-Pin-X-Control Plus

Item no.: BMPXP

The Master Pin X Control Plus system has been specially developed for the fixation of resorbable, non-resorbable, and titanium-reinforced membranes. The pins included are made of grade 5 titanium alloy. The MP30 pins have a self-tapping thread, which allows them to be inserted easily and precisely into the cortical bone. The MP10 pins included have a circumferential groove that increases the surface area of the pins. With their sharp tip and very stable shaft, they can also be inserted easily and precisely into the cortical bone.


INSTRUCTIONS FOR USE

The Bone Management Master Pin X-Control Plus System is indicated for fixing the apical ends of the membrane to the local bone in order to prevent micromobility of the membrane. It enables the fixation of resorbable, non-resorbable, and titanium-reinforced membranes.

APPLICATION

The recommended speeds for use and the maximum permissible speeds are summarized in the product overview under "Contents." These must be observed to ensure proper and safe use.

Setting bleeding points

To refresh the bone, bleeding points can be set in the cortical bone in the region to be augmented using the spiral drills provided. The bone can be pre-drilled using the initial bur if necessary. Pre-drilling is recommended in the case of hard bone.

Application of the augmentation material

In accordance with the respective manufacturer's instructions, the augmentation material is applied to the defect and the selected membrane is adapted to the defect.

Picking up of the MP30 pins

The MP30 pins are located in a sterilizable box that can be opened using a sliding mechanism. The pins are removed from the box using the TL0T1 screwdriver. It is important to ensure a secure grip between the screwdriver and the pin. This is achieved by pressing the screwdriver vertically onto the pin and moving it in an axial rotary motion, first briefly counterclockwise and then clockwise on the pin head. The pins are not fixed in the box. Therefore, care must be taken to avoid vibration when opening the box and picking up the pins with the screwdriver.

Since MP30 pins are implantable products, they must only be picked up using the TL0T1 screwdriver to avoid touching them with your hands.

Picking up the MP10 pins

The pins are stored in a sterilizable box that can be opened using a sliding mechanism. The pins are removed from the box using the pin holder. To do this, press the pin holder vertically onto the pin heads from above. It is important to ensure that the pins are held securely and sit straight in the pin holder. To align the pins straight in the pin holder, the pin holder is ideally held between two fingers when pressing down on the pins and then moved briefly in an axial rotary motion clockwise and counterclockwise on the pin head. This centers the pins in the pin holder. The pins are not fixed in the box. Therefore, care must be taken to avoid vibration when opening the box and when picking up the pins with the pin holder.

Since MP10 pins are implantable products, MP10 pins must only be picked up using the MP12 pin holder to avoid touching them with your hands.

Fixation of the membrane with MP30 pin

To attach the membrane over the defect, the pins are inserted through the membrane into the cortical bone using the TL0T1 screwdriver. The pin can be inserted into the bone without pre-drilling, as it has a self-tapping thread design that prevents aspiration. The pin is gently tightened until the membrane is firmly fixed in place. The bone must not be overloaded. To make it easier to remove the screwdriver from the pin after screwing it in, the screwdriver can be angled at 45°.

Compression of the membrane must be avoided.

Fixation of the membrane with MP10 pin

To attach the membrane over the defect, the pins clamped in the pin holder are applied to the cortical bone through the membrane using a surgical hammer (e.g., M1154), taking care to prevent aspiration. By tilting the pin holder sideways, the holder can finally be released from the applied pin without negatively affecting the fixation of the pin in the bone. Vertical removal of the pin holder should be avoided, as this can exert excessive tensile forces on the pin. In the case of very soft bone, the use of the fixing holder is optionally indicated to fix the pins in the bone when releasing the pin holder. To do this, the uncoated side of the fixing holder is inserted between the lamellae of the pin holder while the pin holder is being detached from the pin.

Compression of the membrane must be avoided.

Removal of the pins

After clinically confirmed successful healing, a small incision is made in the gingiva above the pin to be removed. The MP30 pin can then be carefully loosened using the TL0T1 screwdriver by turning it counterclockwise. The MP10 pin can be carefully loosened using a flat rasp, for example, and removed with tweezers. Alternatively, the MP11 screwdriver can be used for loosening. To do this, turn the pin counterclockwise using the MP11 screwdriver. The pin can then be removed with tweezers. Explantation is also essential if the pins become loose. The pins are a single-use product and therefore cannot be reused. The membrane must be handled in accordance with the manufacturer's instructions.

Sterilization of the pins

The instruments and pins are supplied non-sterile and must be prepared before use (cleaning/disinfection/sterilization). For proper preparation of the items, please follow the instructions for preparation (cleaning, disinfection, and sterilization) of medical devices from Hager & Meisinger GmbH.

Disposal

When disposing of the instruments (at the end of their service life or after the specified shelf life has expired), ensure that the product is disposed of in the waste for biological hazardous substances. All packaging components are disposed of in accordance with national regulations (e.g., dual waste system).

Contraindications

In general, general medical as well as local, absolute, and relative contraindications for dental surgical procedures must be observed.

Absolute contraindications

- Incomplete dentoalveolar growth (exception: cases in which no dentoalveolar growth is to be expected, e.g., ectodermal dysplasia)
- Acute infections in or near the area to be augmented, as well as local pathological processes (e.g., symptoms such as fever, local inflammation, abscesses)
- Insufficient bone quality and/or quantity

Relative contraindications

- Known allergies and/or foreign body reaction to the implant material
- Systemic diseases
- Diseases that impair bone metabolism
- Medical treatments that can impair bone quality and quantity (e.g., cortisone, immunosuppressants)
- Drug and alcohol abuse
- Lack of patient cooperation
- Poor blood circulation
- Heavy labor or active sports
- Mental state that may lead to non-compliance with medical instructions
- Severely atrophied jaw

IMPORTANT: Care must be taken to protect anatomical structures (safety distance of at least 2 mm) and to avoid damage to adjacent teeth/tooth roots (risk of damage, infection/dehiscence).

SPECIAL NOTES FOR PINS: Excessive physical activity and trauma affecting the membrane can lead to premature failure of the pin due to loosening or breakage. The pins may only be used in accordance with the instructions for use and safety instructions. Any modification to the procedure or the product may impair the success of the treatment or damage the product. The pins are a single-use product and must not be reprocessed or reused after use. Even pins that have already been used but are visually undamaged may show microscopic damage due to previous mechanical stress, leading to fatigue and product failure. The pins can be influenced by magnetic fields emanating from an MRI system. Non-clinical tests on equivalent products have shown that titanium pins, are only conditionally MRI-compatible and are therefore only conditionally safe when scanning a patient in an MRI. In these non-clinical tests, using a 3 Tesla MRI system and a gradient-echo sequence, image artefacts were found shadowing an area up to about 2mm around the implant.

Therefore, the following conditions apply when scanning a patient with Pins.

- static magnetic fields of only 1.5 and 3 Tesla
- Gradient fields maximum 3,000 Gauss/cm (30 T/m) (for spatial coding) in normal operation mode:
 - maximum average whole-body MRI system adsorption rate of ≤ 2 W/kg

For a scan duration of 15 min in normal mode under the above conditions, the maximum temperature increase of the pin is <1.5 °C

ASPIRATION HAZARD: The user must ensure that the pins are secured against aspiration both during attachment of the membrane and during removal.

NOTE ON MEMBRANE TECHNOLOGY: Please refer to the information provided by the manufacturer of the membrane you have selected.

Serious incidents related to the product must be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is located.

The Summary of Safety and Clinical Performance (SSCP) provides a summary of clinical data and other information on the safety and clinical performance of the products. The SSCP can be viewed on the website of the European Database for Medical Devices (EUDAMED) <https://ec.europa.eu/tools/eudamed> after publication or can be requested by emailing RA@meisinger.de. Users are encouraged to inform patients about the availability of the SSCP on the aforementioned website.

General information

Please also note the general application and safety instructions for MEISINGER products in the medical field and the instructions for reprocessing (cleaning, disinfection, and sterilization) of medical devices from Hager & Meisinger GmbH at <https://www.meisinger.de/Services/Downloads/>.


						
Manufacturer	Item number	Medical device	Batch number	Do not reuse	Not sterile	Attention
				 ifu.meisinger.de		
Conformity mark	Do not use if packaging is damaged	Date of manufacture	Packaging unit	Follow the electronic instructions for use	Unique product identification	Authorized representative for Switzerland



 **Hager & Meisinger GmbH**
 Hansemannstr. 10
 41468 Neuss • Germany
 Tel. +49 2131 2012-0
 Email info@meisinger.de
 Website: www.meisinger.de

Contents

	Twist drill with stop	Initial bur	RA short screwdriver	Twist drill		Pins	Pins with thread
							
	1	1	1	1	1	20	10
Fig.	203S	186RF	31054	203RF	203RF	MP10*	MP30*
Shaft ¹	204	204	204	205	205	-	-
Size ²	012	018	-	006	008	-	-
Length mm	4.0	12.0	15.5	7.0	7.0	3.65	3.20
Δ mm	-	-	-	-	-	0.87	-
$\square$ mm	1.2	1.8	-	0.6	0.8	0.95	1.45
Opt. speed rpm	1,000	800	10	1,000	1,000	-	-
Max. speed rpm	6,000	1,000	15	2,000	2,000	-	-
MD	CE 0044						

	Screwdriver
	
	1
Fig.	TL0T1
Length mm	148.0
MD	CE

	Pin holder
	
	1
Fig.	MP12
Length mm	170.0
MD	CE

	Screwdriver
	
	1
Fig.	MP11
Length mm	35.0
MD	CE

¹204 = RA 205=RA L ² Largest workpiece diameter in 1/10 mm

  Smallest diameter →  ← Outer diameter


Optionally available

Hammer


	1
Fig.	M1154
Length mm	190.0
	

Screwdriver


	1
Fig.	MP15
Length mm	148.0
	

Fixing holder

	1
Fig.	MP14
Length mm	165.0
	

Pins, flat head

	10
Fig.	MP20*
Length mm	3.15
 mm	0.87
 mm	0.95
Head diameter mm	2.51
	 0044

¹204 = RA ²05=RA L ² Largest workpiece diameter in 1/10 mm

  Smallest diameter →  ← Outer diameter

All previous versions lose their validity with the publication of this instruction manual. New information or safety-related changes compared to the old version are shown in red in the text. Please note that the English original (not the US version) of the instruction manual serves as the template for all translations. In the event of any discrepancies, the English original text is binding.

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