CONTENT

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Application and safety instructions for MEISINGER

myplant bio and myplant two prosthetic components

These instructions for use apply to cover screws, healing abutments, standard abutments, temporary abutments, shoulder abutments, Titanium Bases (Ti-Base), Multi Unit abutments, ball anchors and pre-milled titanium blocks (Pre-Face) of the myplant implant systems. The following information always refers to the two myplant implant systems myplant two and myplant bio. System-specific information is shown separately where applicable.

INTENDED USE

The prosthetic components are indicated for the restoration of dental implants during the healing phase after implantation to enable subgingival and transgingival healing and to enable a prosthetic restoration in the maxilla and mandible for various indications on the basis of a dental implant.

AREA OF APPLICATION

The myplant prosthetic components are suitable for the restoration of myplant dental implants with different lengths and diameters. Prosthetic applications can include customised crowns, bridges, partial or full dentures. Depending on the indication and anatomical conditions, the user utilises the various components for the individual prosthetic purpose.

GENERAL PRODUCT DESCRIPTION

The choice of prosthetic restoration and the prosthetic components required for it have a significant influence on the long-term success of the implantation. Intensive communication between the dentist and dental technician, careful pre-prosthetic planning and patient involvement are important prerequisites for this. The prosthetic restoration can be made using single crowns, bridges, partial or full dentures and can be cemented, bonded, screw-retained or removable on corresponding abutments for myplant implants.

The myplant implant systems offer a wide range of different prosthetic components for the individualised prosthetic restoration of the patient. These include cover screws, healing abutments, standard abutments, temporary abutments, shoulder abutments, Titanium Bases (Ti-Base), Multi Unit abutments, ball anchors and pre-milled titanium blocks (Pre-Face). The connection of the prosthetic components to the implant is ensured via a selflocking, rotation-stable inner cone and, depending on the product variant, via a rotation-secured interface.

PRODUCT-SPECIFIC INSTRUCTIONS FOR THE PROSTHETIC PROCEDURE

The following table shows which prosthetic components are present with or without rotation lock (index).

Prosthetic component	Version with rotation-proof interface (index)	Version without anti-rotation interface
Cover Screw		✓
Healing abutments		✓
Standard abutments	✓	✓
Shoulder abutments		✓
Pre-milled titanium blocks (Pre-Face)	✓	✓
Titanium Bases (Ti-Base)	✓	✓
Multi Unit abutments		✓
Ball anchor		✓

Indexed abutments can only be picked up by an indexed implant interface. Non-indexed abutments can be used with both indexed and non-indexed implant interfaces.

All prosthetic components and instruments are supplied non-sterile and are labelled with the corresponding symbols. They must be cleaned, disinfected and sterilised before use.

All prosthetic components are intended for single-use. If single-use products are reused, there is an increased risk of infection and risk-free functional safety (e.g. in terms of stability or accuracy of fit) cannot be guaranteed.

COMPATIBILITY INFORMATION

Myplant implants must always be used with myplant abutments and prosthetic components from the same implant series (e.g. myplant two implants with myplant two abutments and prosthetic components or myplant bio implants with myplant bio abutments and prosthetic components). For the complete prosthetic restoration, the compatibility of all components and accessories used must always be ensured so that the choice of system components does not impair the restoration (e.g. if an indexed restoration is to be used, only indexed accessories may be used). In addition, care must be taken to ensure that the prosthetic components are only used with the screwdriver provided.

Prosthetic components with and without a rotation-secured interface (index) are available. The non-indexed prosthetics allow free positioning in the implant. The indexed version allows placement in four defined positions and is therefore primarily used to transfer the prosthetic components selected by the dental laboratory from the model back into the mouth in the same position. The special index construction supports self-centring during insertion into the implant. Once the prosthetics have been fixed in the implant using a screw, the internal, self-locking tapered connection ensures that they are secured against rotation and stability.

COVER SCREWS AND HEALING ABUTMENTS

After insertion, the myplant implant is sealed. For subgingival (covered) healing, each implant is supplied with a sterile cover screw with a height of 1.0 mm. Alternatively, a separately available cover screw with a different gingival height (0.0 mm or 2.0 mm) can be used. In the case of subgingival healing, a second procedure is necessary. After the healing phase, the gingiva is surgically opened at the implant site, the cover screw is removed and replaced with a suitable secondary component. For transgingival healing (open), a healing abutment with the desired gingival height should be selected. To shape the soft tissue, the wound edges are sutured atraumatically close around the restoration. The suture closure must not impair the vascularisation of the mucosal edges. The selection of the appropriate closure screw or healing abutment is individualised for each patient and is the responsibility of the user. It is recommended to tighten the cover screw or healing abutment hand-tight (max. 5-7 Ncm).

NOTE: If the cover screw cannot be removed with the screwdriver, the cover screw can be unscrewed from the implant using the unscrewing mandrel (ZAD01/ZAD02) after subgingival healing and exposure.

IMMEDIATE PROVISIONAL RESTORATION

Temporary abutments can be used prior to insertion of the final prosthetic restoration to maintain, stabilise and shape the soft tissue during the healing phase. It is recommended to tighten the temporary abutment with a torque of < 15 Ncm in the implant.

Immediate prosthetic restoration should only be carried out if the patient's prerequisites are met in the early phase of osseointegration, otherwise a delayed solution is preferable. The temporary restoration should not be placed in occlusion and - where possible - the implants should be splinted in this phase.

IMPRESSION

The impression can be taken at implant level using an open or closed impression, at abutment level using an impression tray or digitally, so that a model with laboratory analogues can then be fabricated. When using angled prosthetic components, the position-identical retransfer from the model to the mouth must be ensured (e.g. observing the indexing of the prosthetic components or fabricating a transfer key). When taking impressions at abutment level, prosthetic components that are not secured against rotation must not be removed again after the impression has been taken.

The components used for the impression are screwed hand-tight (max. 5-7 Ncm) into the implant. When using an impression material, care must be taken to select a suitable material. This must be used in accordance with the manufacturer's instructions.

PROSTHETIC RESTORATION

The following table lists the prosthetic treatment options.

Prosthetic component*	Single tooth	Bridge	Partial denture	Total denture
Standard abutment	✓	✓	✓	✓
Titanium base (Ti-Base)	✓	✓		
Multi Unit abutment		✓	✓	✓
Ball anchor			✓	✓
Shoulder abutments	✓	✓		
Pre-milled titanium block (Pre-Face)	✓	✓		

* A single-stage procedure must not be carried out with 6.6 mm implants, as there is insufficient anchoring in the bone. The risk of implant loss is increased for single-tooth restorations with 6.6 mm implants and for free-end bridges. Implants with a length of 6.6 mm are currently only recommended as auxiliary or support implants until corresponding clinical studies are available.

The final prosthetic restoration can be inserted and placed in occlusion as soon as the implant has sufficient stability or is fully osseointegrated.

After receiving the final restoration, the previous temporary restoration is removed. The cleaned and sterilised restoration is then inserted into the patient's mouth. Before tightening the screw, ensure that the implant-abutment connection and alignment are correct.

IMPORTANT: Care should always be taken to ensure even load distribution of the prosthesis.


Torques of more than 15 Ncm (25 Ncm for Multi Unit abutment and ball anchor) can lead to failure of the prosthetic components and/or implant. Torques below the recommended values may lead to loosening of the prosthetic component and thus to failure of the prosthetic component and/or implant. Once the screw has been tightened, it is recommended not to remove it again in order to prevent wear of the screw. If it is necessary to remove the restoration again, it can be reattached with a new screw that matches the prosthetic component.

NOTE: For the correct procedure for bonding or cementing, the manufacturer's instructions of the respective material manufacturer must be observed.

Abutment

The different abutment series enable friction-based, screw-retained, bonded or cemented fixation of the prosthetic components. In the case of angled abutments, in addition to the correct selection of the angulation, the alignment of the angulation must be taken into account.

myplant two superstructures

myplant bio superstructures

Standard abutments

The standard abutments are inserted using the corresponding myplant screwdriver. The indexed abutment can be placed in four different positions. Standard abutments without anti-rotation protection, on the other hand, can be placed freely in the implant. In all cases, care must be taken to ensure that the standard abutment is correctly seated in the implant interface. Taking into account the recommended torque of 15 Ncm, the standard abutment is now tightened in the implant using the insertion tool. It is recommended to check the specified torque again after approximately one week and retighten if necessary. With the myplant standard abutments, screw-retained or cemented prosthetic concepts can be realised in addition to conometric restorations.

Titanium base

The titanium base (Ti-Base) is inserted using the corresponding myplant screwdriver. The indexed abutment can be placed in four different positions. Titanium bases without anti-rotation protection, on the other hand, can be placed freely in the implant. In all cases, care must be taken to ensure that the titanium base is correctly seated in the implant interface. Taking into account the recommended torque of 15 Ncm, the titanium base is now tightened in the implant using the insertion tool. It is recommended to check the specified torque again after about a week and retighten if necessary. The final restoration is cemented onto the residual limb. Cement residues must be thoroughly removed.

For cement-free fixation, the titanium base has an occlusal screw channel through which the prosthesis is screwed occlusally onto the die. The occlusal screw for ceramic is tightened to 10 Ncm using the screwdriver.

Multi Unit abutment

The appropriate Multi Unit abutment is selected according to the implant used and the planned restoration. The Multi Unit abutment is inserted using the corresponding myplant screwdriver. Ensure that the Multi Unit abutment is correctly seated in the implant interface. The Multi Unit abutment is now fixed in the implant using the screwdriver, taking into account the recommended torque of 25 Ncm. Ensure that the prosthesis is seated without tension. It is recommended to check the specified torque again after about a week and retighten if necessary.

As there is no rotation lock between the abutment and superstructure, at least two myplant Multi Unit abutments must always be used together.

Ball anchor

The ball head anchor is inserted using the system-compatible myplant screwdriver. Ensure that it is seated correctly in the implant. The retention ball must be positioned above the mucosa. Taking into account the recommended torque of 25 Ncm, the ball head anchor is now tightened in the implant using the insertion instrument. It is recommended to check the specified torque again after about a week and retighten if necessary.

Shoulder abutments

The shoulder abutment is inserted using the corresponding system-compatible myplant screwdriver. The indexed abutment can be placed in four different positions. Shoulder abutments without anti-rotation protection can be placed freely in the implant. In all cases, care must be taken to ensure that the shoulder abutment is correctly seated in the implant interface. Taking into account the recommended torque of 15 Ncm, the shoulder abutment is now tightened in the implant using the insertion tool. It is recommended to check the specified torque again after about a week and retighten if necessary.

Pre-milled titanium blocks

The pre-milled titanium blocks (Pre-Face) are used to fabricate patient-specific CAD/CAM-milled abutments. All pre-milled titanium blocks are compatible with the Medentika holder and have a rotation-secured implant interface. The appropriate titanium block is selected according to the diameter of the prosthetics to be fabricated. The manufactured abutment is inserted using the corresponding system-compatible myplant screwdriver. The abutment can be placed in four different positions. The alignment of the angulation must be observed. Ensure that the abutment is correctly seated in the implant interface. Taking into account the recommended torque of 15 Ncm, the abutment is now fixed in the implant by tightening the base screw using the screwdriver. It is recommended to check the abutment screw in the implant again after about a week with the specified torque and retighten if necessary.

INTENDED USER PROFILE: Professional users / specialist users in the field of dental implantology.

INTENDED PATIENT POPULATION: Taking contraindications into account, the products can be used for patients of any gender, when the alveolar arch is developmentally mature.

Contraindications

In principle, general medical as well as local, absolute and relative contraindications for dental surgery must be taken into account.

Absolute contraindications include pregnancy (and breastfeeding), serious internal diseases, bleeding disorders, metabolic disorders, metabolic bone disorders (e.g. osteoporosis), reduced immune defence, incomplete growth of the upper and lower jaw and blood clotting disorders. Furthermore, all components are made from materials that are suitable for the stated purposes. In case of allergies or hypersensitivity to the chemical components of the following materials used, the use of the products concerned is contraindicated:

- Titan Grade 4 (3.7065)
- Titan Grade 5 ELI (3.7165 Medical Grade)
- Cerid® 45ON

Relative contraindications include bruxism, psychiatric disorders, medication, drug, alcohol and tobacco consumption as well as insufficient bone supply and poor bone quality, poor oral hygiene and unrestored residual teeth. If the bone supply is insufficient, we refer to the standard augmentation procedures and the relevant literature. Incorrect selection and/or positioning of the prosthetic components can lead to biomechanical overloading due to unfavourable force distribution, resulting in implant loss, fatigue fractures of implants, abutments or prosthetic screws.

Notes

The user must ensure that all components used are secured against aspiration or swallowing.

All prosthetic components are intended for single-use. If single-use products are reused, there is an increased risk of infection and risk-free functional safety (e.g. in terms of stability or accuracy of fit) cannot be guaranteed.

At each step of the prosthetic restoration (e.g. insertion of the cover screw or healing abutment, impression taking and insertion of the final prosthesis), it must be ensured that the internal geometry of the implant is free of contamination so that the implant-abutment connection and thus the success of the restoration are not negatively affected.

All prosthetic components and instruments are supplied non-sterile and are labelled with the corresponding symbols. They must be cleaned, disinfected and sterilised before use. Only system-specific components may be used. Furthermore, any type of modification to the procedure or the product may impair the success of the treatment or damage the product. The patient must be expressly and intensively informed about the necessary oral hygiene measures. Inclusion in a regular recall system is recommended. The myplant prosthetic components have not been evaluated for safety and compatibility in the MRI environment and have not been tested for heating, surface diffusion or image defects in the MRI environment. Safety in the MRI environment is uncertain. Scanning a patient with these products could result in injury to the patient.

Liability and warranty

The use of non-system components and any type of modification may impair the function of the implant system (problems with the fit, reduction in fatigue strength, damage or even fracture of the implant) and excludes any warranty or replacement by the manufacturer. This applies in particular to non-recommended, different application procedures. The processing and application of the products is beyond our control and is the sole responsibility of the user. Any liability for damage caused in this way is excluded. Advice on the technical application of our products is given verbally, in writing, by means of electronic media or by means of demonstrations. It is based on the state of the art in science and technology as known to us at the time the product is placed on the market. It does not release the user from the obligation to personally test the product for its indication-related suitability and to undergo further training in the field of dentoalveolar surgery/implantology.

Serious incidents related to the devices shall be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

Waste disposal

When disposing of the prosthetic components and instruments (at the end of their service life or after expiry of the specified shelf life), ensure that the product is disposed of in the waste for biohazardous substances. All packaging components are disposed of in accordance with national regulations (e.g. dual waste system).


Summary report on safety and clinical performance

The „Summary of Safety and Clinical Performance“ (SSCP) summarises the clinical data and further information on the safety and clinical performance of myplant implants. The SSCP can be viewed on the website of the European Database on Medical Devices (EUDAMED) <https://ec.europa.eu/tools/eudamed> after publication or otherwise requested by e-mail to RA@meisinger.de. The user is required to inform the patient about the availability of the SSCP on the aforementioned website.

General information

Please also observe the general application and safety instructions for MEISINGER products in the medical field and also the Instructions for the reprocessing (cleaning, disinfection and sterilisation) of medical devices from Hager & Meisinger GmbH.

Hager & Meisinger GmbH
Hansestr. 10
41468 Neuss • Germany
Tel.: +49 2131 2012-0
Fax: +49 2131 2012-222
Email: info@meisinger.de
Internet: www.meisinger.de

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